

TRAcelet™

Compression Device



TRANSFORMING YOUR TRANSRADIAL EXPERIENCE

Technical Specifications & Ordering Information

Overview

Working closely with nearly 1000 transradial healthcare providers,¹ we created a comprehensive collection of transradial catheters and accessories to **transform your transradial experience — from access to closure.**

The TRAcelet™ compression device is our innovative solution for simplified pressure control. Its dual-compression balloon system uniformly **delivers pressure to the access site to stop bleeding, while maintaining radial artery patency.**²

Design Rationale

Transforming your transradial experience means finding new approaches to optimise your patient care.

Innovation together — Our collaboration is fuelling purposeful innovation in transradial-specific design. Your input guided our design efforts, and working together, we created an innovative solution for simplified pressure control.

Details matter — The TRAcelet™ compression device integrates key performance characteristics that make a difference — for you and your patients.

A complete solution, completely supported — Our education, training, and support help you optimise your transradial experience.

Features and Benefits

In the cath lab, use the syringe provided with your TRAcelet™ compression device for pressure optimisation.


In the recovery area, use the unique dial-based system for controlled pressure adjustments without the need for a syringe.³

¹ Market research on file at Medtronic.

² Bench test data on file at Medtronic. Bench test data may not be indicative of clinical results.

³ Data on file at Medtronic.

TRAcelet™ Compression Device Technical Specifications



⁴ Data on file at Medtronic.

⁵ Bench test data on file at Medtronic.

Inflation Volume Numbers (cc)	Nominal	13	Strap Lengths (cm)	Regular	19.2
	Maximum	18		Large	25.2

TRAcelet™ Compression Device Ordering Information

CFN	Strap Length
TRACR	Regular
TRACL	Large



Medtronic

Europe
Medtronic Intl. Trading S.A.R.L.
 Route du Molliau 31
 Case Postale
 CH-1131 Tolochenaz
 Switzerland
 Tel: 41.21.802.7000
 Fax: 41.21.802.7900

Canada
Medtronic of Canada Ltd.
 99 Hereford Street
 Brampton, Ontario L6Y 0R3
 Canada
 Tel: 905.460.3800
 Fax: 905.460.3998
 Toll-free: 800.268.5346

Asia Pacific
Medtronic Intl. Ltd.
 49 Changi South Avenue 2
 Nasaco Tech Centre
 Singapore 486056
 Singapore
 Tel: 65.6436.5000
 Fax: 65.6776.6335

Latin America
Medtronic USA, Inc.
 Doral Corporate Center II
 3750 NW 87th Avenue,
 Suite 700
 Miami, FL 33178
 USA
 Tel: 305.500.9328
 Fax: 786.709.4244

medtronic.com